

optipro.TM

surgical masks


the newest addition
to the **OptiPro** range

Optimum Medical's high quality OptiPro Surgical Masks effectively cover the user's nose and mouth in order to minimise the incidence of cross-infection in clinical environments.

Why choose OptiPro Surgical Masks?

- ✓ Our face masks adhere to European Standard EN 14683:2019.
- ✓ Features a protective 3-ply design. This creates an effective barrier, providing superior levels of protection.
- ✓ Stops microbes from entering or exiting the mask, thanks to a melt-blown polymer material which provides a highly effective filter.
- ✓ Comfortable yet secure. The latex-free, elasticated loops slide easily over the ears and securely hold the mask against the user's face. Soft-touch fabric helps to reduce the risk of pressure damage behind the user's ears.
- ✓ The pleat design allows the user to expand the mask and fully cover the area from nose to chin, achieving a perfect fit.

3 LAYERS FOR OPTIMUM PROTECTION



Ordering Information

PR CODE	PRODUCT NAME	INNER PACK	OUTER PACK
1104	OptiPro Surgical Masks	Box of 50	Case of 40

Evidence supports that 3-ply surgical masks are effective in protecting health care workers against transmissible acute respiratory infections in clinical settings. (Smith et al., 2016)

Why choose products from Optimum Medical?

- We're a UK-based manufacturer with years of experience creating high-quality medical devices.
- Our products are used in every NHS hospital.
- We export our products to over 60 countries around the world.
- We're regulated to ISO 13485 international standard.

To find out more information, email:
enquiries@optimummedical.co.uk

Or give us a call directly on:
+44 (0) 113 263 3849

References:
Smith, J. D., MacDougall, C. C., Johnstone, J., & Copes, R. A. (2016, May 17). Effectiveness of N95 respirators versus surgical masks in protecting health care workers from acute respiratory infection: a systematic review and meta-analysis. *Canadian Medical Association Journal*, 188(8), 567-574. Retrieved April 24, 2020



optimummedical® | constantly creating better

Tennant Hall, Blenheim Grove, Leeds, LS2 9ET, UK